

FIX Voxel size: 0.9×0.9×0.9 mmPAT: 3 Rel. SNR: 1.00 : nih5i

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Position	R1.4 A22.9 H20.6 mm
Orientation	T > C-14.1
Phase enc. dir.	A >> P
AutoAlign	---
Slab Scale	-10 %
Slices per slab	12
FoV read	190 mm
FoV phase	100.0 %
Slice thickness	0.90 mm
TR 1	63.5 ms
TR 2	8867 ms
TE 1	8.75 ms
TE 2	28.55 ms
TE 3	48.35 ms
Averages	1
Filter	None
Coil elements	A32

Contrast - Common

TR 1	63.5 ms
TR 2	8867 ms
TE 1	8.75 ms
TE 2	28.55 ms
TE 3	48.35 ms
Multi-echo spacing	19.80 ms
Magn. preparation	Non-sel. HSN IR
TI 1	713 ms
TI 2	1475 ms
TI 3	2237 ms
Flip angle	45 deg
Fat suppr.	Fat sat.
Magn. Prep. Shots	3

Contrast - Dynamic

Averages	1
Averaging mode	Long term
Reconstruction	Magn./Phase
Measurements	80
Pause after meas.	0.0 s

Resolution - Common

FoV read	190 mm
FoV phase	100.0 %
Slice thickness	0.90 mm

Resolution - Common

Base resolution	212
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	6/8
Slice partial Fourier	Off
Interpolation	Off

Resolution - iPAT

PAT mode	CAIPIRINHA
Acc. factor PE	1
Ref. lines PE	80
Acc. factor 3D	3
Ref. lines 3D	12
CAIPI 3D Shift	1
Reference Scan Mode	EPI/separate
CAIPI Mode (tooltip)	Skipped-CAIPI
Total PAT factor	3

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slab group	1
Slabs	1
Position	R1.4 A22.9 H20.6 mm
Orientation	T > C-14.1
Phase enc. dir.	A >> P
Slab Scale	-10 %
Slices per slab	12
FoV read	190 mm
FoV phase	100.0 %
Slice thickness	0.90 mm
TR 1	63.5 ms
TR 2	8867 ms

Geometry - AutoAlign

Slab group	1
Position	R1.4 A22.9 H20.6 mm
Orientation	T > C-14.1
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R1.4 A22.9 H20.6
R	1.4 mm
A	22.9 mm
H	20.6 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-14.1
> S	0.0

Geometry - Saturation

Saturation mode	Standard
-----------------	----------

SIEMENS MAGNETOM Investigational_Device_7T_Plus

Geometry - Saturation

Fat suppr.	Fat sat.
------------	----------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	R3.8 A20.2 H19.9 mm
! Orientation	T > C-14.2
! Rotation	0.35 deg
! A >> P	211 mm
! R >> L	117 mm
! F >> H	40 mm
Reset	Off

System - Tx/Rx

Frequency 1H	297.204286 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
! Ref. amplitude 1H	170.000 V

Sequence - Part 1

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	3
Echo spacing	1.01 ms
Bandwidth	1180 Hz/Px

Sequence - Part 2

EPI factor	18
Segmentation	9
RF pulse type	Normal
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On
Turbo factor	12

Sequence - Special

PATRef FA	3 deg
RF duration	1200 us
RF BWT product	12
Ernst T1	1200 ms
PATRef prep. shots	10
Volume dummy shots	0
Dummy Measurements	0
ETL per RTEB	1
Invert PE	Off
Min. TE if PF	On
Echo Time Shift	On
Ramp Sampling	On
NORDIC	Off
SVDPC	On
Sym VASO	Off
Dual-pol. EPI	Off
Invert RO	Off
Invert 3D	Off
Disable PF reco	Off
Disable PF reco	Off
Save sampling	Off
PE VComp	Off
Water Exc.	-none-
External PC	per Series
Saturation RF	per Shot
FIDNavs	-none-
EPI rise time factor	1.10
Mosaic DICOMs	On
Modify Ice Config	On
G-factor map	Off
GRAPPA Regularization	50000 10^-6
HSN RF power scale	3.00
Inversion Delay	322 ms
Relaxation Delay	0 ms
Var. FA /MAGEC	0

Sequence - Assistant

Mode	Off
------	-----

SIEMENS MAGNETOM Investigational_Device_7T_Plus

Geometry - Saturation

Fat suppr.	Fat sat.
------------	----------

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Brain
B1 Shim mode	TrueForm
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	R3.8 A20.2 H19.9 mm
! Orientation	T > C-14.2
! Rotation	0.35 deg
! A >> P	211 mm
! R >> L	117 mm
! F >> H	40 mm
Reset	Off

System - Tx/Rx

Frequency 1H	297.204286 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
! Ref. amplitude 1H	170.000 V

Sequence - Part 1

Introduction	On
Dimension	3D
Reordering	Linear
Contrasts	3
Echo spacing	1.01 ms
Bandwidth	1180 Hz/Px

Sequence - Part 2

EPI factor	18
Segmentation	9
RF pulse type	Normal
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On
Turbo factor	12

Sequence - Special

PATRef FA	3 deg
RF duration	1200 us
RF BWT product	12
Ernst T1	1200 ms
PATRef prep. shots	10
Volume dummy shots	0
Dummy Measurements	0
ETL per RTEB	1
Invert PE	Off
Min. TE if PF	On
Echo Time Shift	On
Ramp Sampling	On
NORDIC	Off
SVDPC	On
Sym VASO	Off
Dual-pol. EPI	Off
Invert RO	Off
Invert 3D	Off
Disable PF reco	Off
Disable PF reco	Off
Save sampling	Off
PE VComp	Off
Water Exc.	-none-
External PC	per Series
Saturation RF	per Shot
FIDNavs	-none-
EPI rise time factor	1.10
Mosaic DICOMs	On
Modify Ice Config	On
G-factor map	Off
GRAPPA Regularization	50000 10^-6
HSN RF power scale	3.00
Inversion Delay	322 ms
Relaxation Delay	0 ms
Var. FA /MAGEC	0

Sequence - Assistant

Mode	Off
------	-----